

# THORNE

R E S E A R C H



## KRILL OIL

**60 FISH GELATIN GELCAPS**

**DIETARY SUPPLEMENT**

## Supplement Facts

Serving Size: Two Gelcaps  
Servings Per Container: 30

| Two Gelcaps Contain:    |        | % DV | Two Gelcaps Contain:                                                                          |        | % DV |
|-------------------------|--------|------|-----------------------------------------------------------------------------------------------|--------|------|
| Calories                |        | 8.6  | DHA                                                                                           |        |      |
| Calories from Fat       |        | 6.4  | (Docosahexaenoic Acid)                                                                        | 60 mg  | **   |
| Total Fat               | 0.8 g  | 1%*  | Phospholipids                                                                                 | 400 mg | **   |
| Cholesterol             | 16 mg  | 6%*  | *Percent Daily Values are based on 2,000 calorie diet.<br>**Daily Value (DV) not established. |        |      |
| EPA                     |        |      |                                                                                               |        |      |
| (Eicosapentaenoic Acid) | 150 mg | **   |                                                                                               |        |      |

Other Ingredients: Gelatin (from tilapia), Purified Water and Glycerin (vegetable source) gelcap, Vitamin E (mixed tocopherols).

Contains ingredients derived from shellfish (krill) and fish (gelatin = tilapia).

Suggested Use: Take 1-2 gelcaps two times daily or as recommended by your health-care practitioner.

Tamper Evident: Use only if bottle is sealed.

If pregnant, consult your health-care practitioner before using this product.

Store tightly sealed in a cool, dry place.



Manufactured For: Thorne Research, Inc.  
P.O. Box 25, Dover, Idaho 83825 USA  
1-800-228-1966 [www.thorne.com](http://www.thorne.com)



SP612  
LSP61206