RECOMMENDED USE: Take 2 cansules on workout days. During intense dieting take 1 capsule on non-workout days before bedtime. Consume CortSTAT with at least 6-8 nunces of water

STORAGE INSTRUCTIONS: Protect from heat and freezing. Keep away from local snare heat sources such as snare heaters and radiators. Do not evenue to direct sunlight or leave in vehicle exposed to sunlight. Do not allow exposure to high humidity (above 65%). Keep container tightly closed when not in use.

DESCRIPTION: Flexated cortisol levels in the athlete are known to break down. muscle tissue via inhibiting muscle protein synthesis. CortSTAT by SAN is a bona fide new formula that supports the maintenance of muscle tissue nost-workout." Very recent research has shown that lowered levels of cortisol enable creater amino acid utilization for enhanced nitrogen retention, glycogen synthesis and protein synthesis.* As a dietary supplement, use CortSTAT as a key substrate to

thwart exercise induced stress muscle tissue cannibalization

prevent any disease.

SAN products are manufactured under GMP (Good Manufacturing Practices) standards with stringent quality control quidelines. Additionally each lot produced by SAN is independently tested for purity and potency by HPLC analysis to ensure

only the highest quality ingredients in our innovative formulas."

Manufactured in a facility that also processes egg, shellfish, fish, soy, peanuts, tree nuts and milk ingredients * These statements have not been evaluated by the Food and Drug Administration. This product is not intended to diagnose, treat, cure or SAN* Live a VIBRANT LIFE!



POST-WORKOUT TRAINING **MODULATOR***











Supplement Facts Serving Size: 2 Capsules

200 mg**

Servings Per Container: 30

Alanyl-L-Glutamine Dipeptide (GlutaDiPEP") 800 mg**

Amount Per Serving

Phosphatidylserine 50% (SerinAid®1

Daily Value not established SerinAid® is a registered trademark of Chemi Nutra LLC.



