## **Supplement Facts**

Serving Size: 1 scoop (12 grams)

Serving Per Container: 25-50

Amount Per Serving % Daily Value

Proprietary Blend with Micro-Rx™ & Extend-Rx™ Technologies: 12-24 grams Disodium Creatine Phosphate Tetrahydrate, Creatine Phosphate (CreaPhos™ brand), Micronized Creatine Monohydrate, Creatinol-O-Phosphate, Beta-Alanine Ethyl Ester, L-Taurine Ethyl Ester, Creatine Ethyl Ester, Dicreatine Malate, Inosine Ethyl Ester Adenosine Triphosphate Sodium, Guanidinopropionic Acid, Creatine Alpha-Ketoglutarate, Stearic Acid

\*Daily Value Not Established

Other Ingredients: Naturally Sweet™ Brand Sweetener, Erythritol, Maltitol, Maltodextrin, Xylitol, Tagatose, Sucralose, Silica, Natural and Artificial Flavoring

Suggested Use: Take 1-2 scoops of Anavar®-PCr 30 minutes prior to workout to increase Creatine and ATP saturation with 8-10 ounces of water. Should be taken even on non-training days.



HI-TECH PHARMACEUTICALS, INC. 6015-B Unity Drive • Norcross, GA 30071 www.hitechpharma.com • 1.888.855.7919

Anavar®-PCr is the world's first and only product to feature Disodium Creatine Phosphate Tetrahydrate, which is a proprietary Creatine Series Analog consisting of a novel complex sodium salt of Creatine Phosphate that enhances water solubility. Once in the bloodstream, a sodium chloride dependant transporter is responsible for carrying the creatine molecule directly to muscle cells for ATP synthesis. The bonding of Creatine with Disodium Phosphate Tetrahydrate is the most significant breakthrough in the uptake and absorption of Creatine. † Scientists have known for years that Creatine Phosphate (PCr) elevation is the key for increased energy and fuel to muscle, increased hydration, increased muscular endurance and workload capability, and with decreased muscle fatigue because of accelerated replenishment of cellular Adenosine Triphosphate (ATP). In short, Anavar®-PCr is a revolutionary metabolic catalyst that drastically amplifies Creatine Phosphate (PCr) levels unlike any other product- which leads to explosive increases in muscle mass and strength!+

Anavar®-PCr utilizes Micro-Rx™ and Extend-Rx™ technologies to enhance bioavailability and deliver creatine in "Sustained-Release" microparticles. These microparticles release the product gradually in a consistent manner over a prolonged period of time. Hi-Tech's proprietary process consists of heating and mixing USP Stearic Acid, which is a mixture of Stearic and Palmitic Acids with the active components of Anavar®-PCr until well coated, and then cooling to room temperature until hard. The resultant particles are small, free-flowing, and exhibit release profiles that are released over 12 hours! Anavar®-PCr won't cause the stomach distress associated with other creatine products due to creatine degradation. Anavar®-PCr will render all other creatine products obsolete. Anavar®-PCr is the most effective performance-enhancing muscle-building supplement ever developed!

> †These statements have not been evaluated by the Food and Drug Administration. This product is not intended to treat, cure, diagnose or prevent any disease.

## NAVAR-P(F

## World's Strongest PhosphoCreatine and ATP Elevator

- · First Disodium Creatine Phosphate Tetrahydrate product which delivers Creatine directly into the muscles: • The only ATP & Creatine Phosphate product that is Based Upon Scientifically-Engineered
  - "Sustained Release" for all day Creatine Saturation: The Biggest Breakthrough in Sports Supplementation and Scientifically-Formulated for Powerful Anabolic Muscle Growth and Explosive Power!+

Revolutionary Muscle-Building Catalyst Designed to Rapidly Increase Strength & Size

Tropical Punch Flavor