

Re-Order No. 700149

MAJOR[®] Poly-Vita Drops™ with Iron

COMPARE TO
VITAMINS IN
Poly•Vi•Sol[®]
with Iron

MULTIVITAMIN SUPPLEMENT FOR INFANTS AND TODDLERS

- Delicious Natural Fruit Flavor
- Alcohol and Sugar Free

1 2/3 fl. oz. (50 mL)

USE ONLY WITH ENCLOSED SYRINGE

DIRECTIONS: Administer this supplement using only the device provided. Shake well and fill device to the 1.0 mL line. Dosage not to exceed 1.0 mL per day. Dispense directly into mouth or mixed with formula, fruit juice, cereal or other food. Rinse and retain device for future dosing of this product. **USE FULL DOSAGE.**

Supplement Facts

Serving Size 1.0 mL (Syringe)

Amount per mL	% Daily Value	
	Infants under 12 months	Children under 4 years
Vitamin A (as palmitate)	1500 IU	100%
Vitamin C (as ascorbic acid)	35 mg	100%
Vitamin D (as cholecalciferol)	400 IU	100%
Vitamin E (as dl-alpha tocopheryl acetate)	5 IU	100%
Vitamin B1 (as thiamine HCl)	0.5 mg	100%
Vitamin B2 (as riboflavin-5-phosphate)	0.6 mg	100%
Niacin (as niacinamide)	8 mg	100%
Vitamin B-6 (as pyridoxine HCl)	0.4 mg	100%
Iron (as ferrous sulfate)	10 mg	67%
		100%

Other ingredients: Purified Water, Polysorbate 80, Propylene Glycol, Glycerin, Natural Cherry Flavor, Caramel Color, Potassium Sorbate (Preservative), Sodium Benzoate (Preservative), Citric Acid or Sodium Hydroxide may be used to adjust pH.

WARNINGS: While taking this product, a slight darkening of the teeth may occur. Brushing will minimize this temporary condition. Accidental overdose of iron-containing products is a leading cause of fatal poisoning in children under six. **KEEP THIS PRODUCT OUT OF REACH OF CHILDREN. In case of accidental overdose, call a doctor or Poison Control Center immediately.**

Store at room temperature 15°-30°C (59°-86°F). Protect from freezing. After opening, store away from direct light.

**TAMPER EVIDENT: DO NOT USE IF SAFETY SEAL
UNDER CAP IS BROKEN OR MISSING**

Questions or comments? Call (800) 616-2471

Distributed by: **MAJOR[®] PHARMACEUTICALS**

31778 Enterprise Drive, Livonia, MI 48150 USA
M-133 Re-Order No. 700149

CT04050L Rev. 02/13