

Suggested Use: To use **Liquid Zinc Assay** as a supplement, take 1 teaspoon in water or juice daily or as directed by a health professional.

To use as a zinc test: hold 1 teaspoon in the mouth for 15 seconds. If an immediate perception of taste occurs, the person may have an adequate zinc status. If a lack of taste or delay in taste perception occurs, it may indicate a zinc deficiency. **Liquid Zinc Assay** may be swallowed after testing.

► Our liquid zinc is made using 100% purified water, not filtered or tap water which may contain impurities.



Keep out of reach of children.

****THESE STATEMENTS HAVE NOT BEEN EVALUATED BY THE FOOD AND DRUG ADMINISTRATION. THIS PRODUCT IS NOT INTENDED TO DIAGNOSE, TREAT, CURE OR PREVENT ANY DISEASE.**

R18-0521 R5.3



Liquid Zinc Assay

Highly Absorbable Liquid Zinc**

Rapid Zinc Status Analysis

A DIETARY SUPPLEMENT

8 FL OZ (235 mL)

- ✓ PRL's Quality Guarantee
- ✓ Pure Vegan
- ✓ NON-GMO
- ✓ Purity Verified
- ✓ Potency Verified
- ✓ Tested for Added Adulterants
- ✓ No Chemical or Radiation Sterilization
- ✓ Pathogen Microbiology Performed
- ✓ Heavy Metal Tested
- ✓ Identity by Infra-red Spectroscopy

Supplement Facts

Serving Size 1 Teaspoon (5 mL)
Servings Per Container about 47

Amount Per Serving	%Daily Value
Zinc (Zinc Sulfate Heptahydrate)	10 mg 90%

Other Ingredients: Purified Water

PRL's Quality Program. Manufactured to avoid: added magnesium stearate (an undesirable excipient), artificial colors, artificial flavors or artificial preservatives.

Manufactured in a cGMP compliant facility for superior quality assurance

Tamper Seal: Use only if seal is intact. Consult your health care practitioner if you are nursing/lactating, taking medications or have a medical condition, before taking this or any other product. Store in a cool, dry place.