

Duosoluble™

QGel® Ultra

COENZYME Q₁₀ 60 mg

+ 100 mg Vitamin E

**Cellular Energizer
Enhanced Absorption†**

Dietary Supplement

**500
Softgels**

Directions: Adults: Take one (1) softgel daily as a dietary supplement or as directed by your healthcare professional. Softgel to be taken with or after a

Supplement Facts

Serving Size: 1 Softgel

	Amount Per Softgel	% Daily Value*
Calories	6	
Vitamin E (as d-alpha tocopheryl acetate USP/NF)	100 mg	667%
Coenzyme Q ₁₀ (ubidecarenone USP)	60 mg	**

* Percentage of Daily Value for Adults and Children 4 or more years of age.

** Daily Value not established

Other Ingredients: Polysorbate 80, gelatin, triacetin, medium chain triglycerides, glycerin, sorbitol, high PC lecithin, soybean oil, purified water, mannitol, titanium dioxide and annatto extract.

Contains soy

Manufactured and packaged in NY & MD USA with USA and imported ingredients. By: **Gel-Tec**, Division of **Tishcon Corp.**

30 New York Ave, Westbury, NY 11590

www.CoQ.com

Free of: milk or milk by-products, egg or egg by-products, fish or fish by-products, shellfish or shellfish by-products, tree nuts, wheat or wheat by-products and peanuts or peanut by-products.

KEEP OUT OF THE REACH OF CHILDREN. If you are pregnant or nursing a baby, seek the advice of a health professional before using this product.


THE BOTTLE CAP ON THIS PRODUCT IS SECURED WITH AN IMPRINTED PLASTIC SEAL. DO NOT USE THIS PRODUCT IF SEAL IS CUT, TORN, BROKEN OR MISSING.

Keep bottle tightly closed. Store between 15°-30°C (59°-86°F).

Do not refrigerate

†This statement has not been evaluated by the Food and Drug Administration. This product is not intended to diagnose, treat, cure or prevent any disease.

Manufactured by the BioSolv® process. BioSolv®, Q-Gel® and CoQPure logo are registered trademarks of Tishcon Corp. Duosoluble™, Hydrosoluble™ and Liposoluble™ are trademarks of Tishcon Corp. To report a serious adverse event or to obtain product information, contact 1-800-848-8442.



Control No.:
Exp. Date

7660-SEG
R1120