

VITAL-ZYMES CHEWABLE

COMPREHENSIVE ENZYME SUPPLEMENT

> 180 CHEWABLE **TABLETS**

Supplement Facts

Serving Size 2 Tablets - Servings Per Container 90

| mount Per 2 Tablets | % Daily Value | Amount Per 2 Tablets | % Daily Val |
|--|---------------|--|-------------|
| omprehensive Enzyme Blend | 320 mg* | Phytase | 4 PUΔ |
| With the following enzyme activity: | | Beta-Glucanase | 8 BGUA |
| Carbohydrate Specific Enzymes | | Alpha-Galactosidase | 40 AGSUA |
| Amylase | 732 BAUA | Galactomannase | 80 HCUA |
| Alpha-Amylase | 40.DUΔ | Protein and Peptide Specific Enzymes | |
| Glucoamylase (with isomaltase side chain activity) | 12 AGUA | Acid Protease | 20 SAPUA |
| Sugar Specific Enzymes | | Alkaline Protease | 20 PCA |
| Lactase | 600 ALUA | Protease/Peptidase Complex with endo- | 16,000 HUTA |
| Sucrase (invertase) | 160 SUA | endopeptidase, exopeptidase, and DPP-IV activity | |
| Maltase (diastase) | 80 DP*∆ | Fat Specific Enzyme | |
| Pullulanase (debranching enzyme with isomaltase side chain activity) | 10 UA | Lipase | 600 FIPA |
| | | Other Enzymes | |
| Vegetable/Plant Fiber Specific Enzymes | | Serratia peptidase (enteric-coated): | 1,000 SPU |
| Cellulase | 200 CUA | March 1997 | |
| Hemicellulase/Pectinase/Phytase Complex | 80 HSUA | Daily Value not established. | |

Other ingredients: Xylltol, mannitol, stearic acid, silicon dioxide, and natural cherry flavor.

Manufactured by: SFI USA • Reno, NV 89521 • Klaire.com Caution: If you are pregnant, nursing, have a medical condition, or taking prescription drugs, consult your healthcare professional before using this

product. Keep out of reach of children. Free of the following common allergens: milk/casein, eggs, fish, shellfish, tree nuts, peanuts, wheat, gluten, and soybeans. Contains artificial colors, flavors, or preservatives.

Suggested Use: 2 tablets at the beginning of each meal or as directed by a healthcare professional.

AUnits are those used in edition nine (IX) of the Food Chemicals Code: FCC Peptizyme SP*, a registered trademark of Specialty Enzymes. Do not use if shrinkwrap is broken or missing. Store in a cool, dry place

(59°F-85°F) away from direct light. Produced under a strict quality management 2 rev17072451 system in compliance with Good Manufacturing Practices (GMPs) and

third-party quality certifications.