

SUGGESTED USAGE: As a dietary supplement, take 1 soft-gel every 3 days with a meal containing fat, or as directed by your healthcare practitioner.

HOW IT WORKS: The major biologic function of Vitamin D is to maintain healthy blood levels of calcium and phosphorus. Vitamin D supports bone health by promoting calcium absorption and bone mineralization.* Vitamin D also plays an important role in immune function and helps to regulate cell growth and differentiation.*

CAUTIONS/INTERACTIONS: Although the Upper Limit for Vitamin D is 4,000 IU per day, comprehensive research demonstrates that higher dosages are safe. However, this product should be taken under the supervision of a healthcare professional. Vitamin D should be used cautiously by those taking cardiac glycosides and anyone with hypercalcemia. Vitamin D supplements may interact with atorvastatin. If you are taking any of these prescription medications, consult your healthcare practitioner before using this product. For more information on this product go to www.protocolforlife.com/P0376

*These statements have not been evaluated by the Food and Drug Administration. This product is not intended to diagnose, treat, cure or prevent any disease.

PROTOCOL
FOR LIFE BALANCE®

High Potency Vitamin D₃ 10,000 IU

- Helps Maintain Strong Bones*
- Supports Immune System*

120 Softgels

A Dietary Supplement



Supplement Facts

Serving Size 1 Softgel

	Amount Per Serving	% Daily Value
Vitamin D ₃ (as Cholecalciferol) (from Lanolin)	10,000 IU	2500%

Other ingredients: Extra Virgin Olive Oil and Softgel Capsule (gelatin, glycerin, water).

Formulated & Distributed by Protocol For Life Balance®
395 S. Glen Elyn Rd., Bloomingdale, IL 60108, USA

Not manufactured with yeast, wheat, gluten, milk, egg, fish, shellfish or tree nut ingredients. Produced in a GMP facility that processes other ingredients containing these allergens.

This bottle contains an AGELESS® Oxygen Absorber to ensure freshness. Do not eat packet. Store in a cool, dry place after opening. Please Recycle.

CODE P0376 V1

