As a dietary supplement, mix 1-2 scoops of Gorilla Mode Nitric in 12-14 ounces of water and consume 30 minutes place in 2-14 ounces of water and consume 30 minutes prior braining. Vary the amount of water to achieve your desired flavor level. First time users should begin use with 1/21 scoop or less to evaluate tolerance. DO NOT EVCED 2 SCOOPS IN ANY 24 HOUR PERIOD.

# KEEP OUT OF REACH OF CHILDREN.

This product should not be taken by individuals who are pregnant, nursing, have or suspect a medical condition. are taking any medications or are under 18 years of age. Consult a physician before taking this or any other dietary supplement. Do not use this product with PDE5 inhibitors, nitrites, or other blood pressure lowering medications and/or supplements. Store in a cool, dry

Fiftese statements have not been evaluated by the Food and Drug Administration. This product is not intended to diagnose, treat, cure, or prevent any disease.











MODE NITRIC

STIMULANT FREE PRE-WORKOUT FORMULA



RAISES NITRIC OXIDE: • INTRAMUSCULAR HYPER-HYDRATION: • INCREASES POWER & ENDURANCE



# Supplement Facts

Serving Size: 1 Scoop (17.375g)


1 Scoop (17.375g) Amount Per Serving %DV			2 Scoops (34.75g) Amount Per Serving % (I)		
Sodium (from sodium nitrate) Potassium (from inositol-stabilized arg	202.5 mg 7 37.5 mg 1	%'	Sodium (from sodium nitrate) Potassium (from inositol-stabilized arg	405 mg 1	49
L-Citrulline Creatine Monohydrate Betaine Anhydrous	5000 mg 2500 mg 2000 mg	† † †	L-Citrulline Creatine Monohydrate Betaine Anhydrous	10000 mg 5000 mg	1 1
GlycerPump™ (65% Glycerol Powder)	2000 mg	t	GlycerPump™ (65% Glycerol Powder)	4000 mg 4000 mg	
Malic Acid Agmatine Sulfate	1500 mg 750 mg	†	Matic Acid Agmatine Sulfate	3000 mg 1500 mg	
Nitrosigine® 750 mg † (inositol-stabilized arginine silicate)		Nitrosigine® (inositol-stabilized arginine s	1500 too		
Sodium Nitrate	750 mg	†	Sodium Nitrate	1500 mg	
VasoDrive-AP® 127 mg †   soleucyl-prolyl-proline (IPP) and Valyl-prolyl-Proline (VPP) isolated from hydrolyzed milk casein)		VasoDrive-AP* (Isoleucyl-prolyl-proline (IPP) and Valyl-prolyl-Proline (VPP) isolated from hydrolyzed milk caseio)			

Note: Due to the very high concentration of

stured Exclusively For: Gorilla Mind, LLC 7154 W State Street, STF 138, Roise, ID 83714, U.S.A.

GLYCER PUMP

GlycerPump™ is a trademark of Pinnacle Ingredients (Australia)





Nitrosigine® and the Nitrosigine® logo are registered trademarks of Nutrition 21, LLC.

