

SUGGESTED USE: Take 1-3 capsules two times daily or as recommended by your health-care practitioner.

Tamper Evident: Use only if bottle is sealed.

Store tightly sealed in a cool, dry place.

If pregnant, consult your health-care practitioner before using this product.

Manufactured for:

Tesseract Medical Research, LLC

1370 Trancas St., #376, Napa, CA 94558

www.TessMed.com

Product: GI000

ProButyrate®

Sustained Release
Hypoallergenic



Tesseract
Medical Research®

Dietary Supplement
120 Vegetarian Capsules

SUPPLEMENT FACTS

Serving Size: One Capsule

Servings Per Container: 120

Each Capsule Contains:		% DV
Butyric Acid Complex**	300 mg	*

*Daily Value (DV) not established.

Other Ingredients: CyLoc™ alpha dextrin fiber matrix Butyric acid complex 300 mg.*, high and low viscosity hydroxypropyl methylcellulose (sustained release), hypromellose (derived from cellulose) capsule, calcium laurate, leucine, and plant amylase. *Stabilized Butyric acid as soluble, pure, individual free fatty acid in a hexadextrin matrix molecular trap compounded with DexKey™ release agent (patent pending, Al Czap, Czap R&D, licensed by Tesseract Medical Research, LLC). CyLoc™ and DexKey™ are trademarks of Tesseract Medical Research, LLC.