

SUGGESTED USE: Take 1 capsule two times daily or as recommended by your health-care practitioner. Tamper Evident: Use only if bottle is sealed. Store tightly sealed in a cool, dry place. If pregnant, consult your health-care practitioner before using this product.

Manufactured for:
Tesseract Medical Research, LLC
1370 Trancas St., #376, Napa, CA 94558
www.TessMed.com

RLA Pro SR™

R-Lipoic Acid Enhanced Delivery
Sustained Release
Hypoallergenic



Tesseract
Medical Research®

Dietary Supplement
60 Vegetarian Capsules

SUPPLEMENT FACTS

Serving Size: One Capsule

Servings Per Container: 60

Each Capsule Contains:	% DV
R-Lipoic Acid Proprietary Complex†	250 mg*

*Daily Value (DV) not established.

Other Ingredients: Hypromellose (derived from cellulose) capsule, high and low viscosity hydroxypropyl methylcellulose (sustained release), calcium laurate, and plant derived enzyme.

†Solubility enhanced R-Lipoic Acid in CyLoc™, a dextrin fiber matrix molecular trap. CyLoc™ complexes compounded with DexKey™ release agents (patent pending, Al Czap, Czap R&D, Licensed by Tesseract Medical Research, LLC).