

Vitargo® S2 is the original super carb, proven in university studies in humans to be up to 2.3x faster than maltodextrin (homopolysaccharide). Vitargo's patented IVg technology delivers faster gut transit, glycemic and insulin responses, muscle glycogen refueling, and performance. No other carb has this span of proof. Vitargo S2 is 100% sugar-free, lab tested gluten-free, and certified undetectable banned substances.

- **Leaves the stomach quicker** 130% faster than maltodextrin + sugars – in the first 10 minutes after ingestion.¹ This leads to less stomach "distress" and faster absorbing muscle energy.
- **Gets into muscle faster** 68% faster glycogen re/fueling than maltodextrin + sugars after intense, exhaustive workouts.²
- **Boosts performance in your next workout** Up to 23% greater maximal endurance (average of 10% greater) 2 hours after exhaustive, glycogen-depleting exercise, compared to maltodextrin + sugars.³
- **Turns off muscle protein breakdown** 178% greater insulin than maltodextrin + sugars, within 10 minutes³ – the most potent, natural way to activate the anti-catabolic signals that spare muscle protein.


GENR8-Specific Research — this actual product is university proven in HUMANS.

Proof Before Promises®

References

1. Leiper JB, et al. Improved gastric emptying rate in humans of a unique glucose polymer with gel-forming properties. *Scand J Gastroenterol* 2000; 35:1143-9.
2. Aulin KP, et al. Muscle glycogen resynthesis rate in humans after supplementation of drinks containing carbohydrates with low and high molecular masses. *Eur J Appl Physiol* 2000; 81:346-51.
3. Stephens FB, et al. Post-exercise ingestion of a unique, high molecular weight glucose polymer solution improves performance during a subsequent bout of cycling exercise. *J Sports Sci* 2008; 26:149-54.

These statements have not been evaluated by the Food and Drug Administration. This product is not intended to diagnose, treat, cure, or prevent any disease.



2X Faster Glycemic Rise
1.7X Faster Glycogen

PRE • INTRA • POST
VITARGO® S2™

Fastest Muscle Fuel™

- Faster Absorbing Muscle Energy
- Up to 23% Greater Post-Workout Performance
- 178% Greater Insulin Release for Anti-Catabolic Action

Natural Juicy Orange Flavor

SUGAR-FREE DRINK MIX *Not a low calorie product*

Dietary Supplement Net Wt. 1.7 lbs (770g)

Gluten-Free



PATENTED IVg TECHNOLOGY™

INSTRUCTIONS

Vitargo® S2™ is a different engineered carbohydrate – please follow directions closely:

HOW In a shaker bottle

- Fill a large shaker bottle half full (10 - 12 oz.) with water (room temp. is best).
- Add 2 level scoops of Vitargo® S2 and shake vigorously for 10-15 seconds.
- Add cold water to top and shake vigorously for 15 - 20 seconds.

In a blender

- Vitargo S2 can be combined with your favorite post-training blender concoction – protein powder, branched chain amino acids, etc. and your favorite liquid.
- For best results add 2 level scoops to a blender with water while it's on and blending.

WHEN Before, between, or during exercise/competition

- Mix 1 or 2 level scoops of Vitargo S2 as described above. Note that Vitargo S2's unique molecular viscosity makes it thicker than other carbohydrate drinks – experiment and see if you prefer 1 or 2 scoops. Drink to rapid fueling!

Post-training/competition/glycogen loading

- Take 2 level scoops as soon as possible after finishing. For optimal recovery, take another 1 - 2 scoops 30 - 90 minutes after your first serving. For maximal glycogen loading, take 3 servings/day.

For more detailed info, go to www.genr8speed.com. IF YOU DON'T TRAIN, DON'T TAKE THIS PRODUCT.

Supplement Facts

Serving Size: 2 level scoops
Servings Per Container: 10

	Amount Per Serving	%DV*
Calories	280	
Calories from Fat	0	
Total Fat	0 g	0%
Total Carbohydrate	70 g	23%
Sugars	0 g	
Protein	0 g	0%

*Percent Daily Values (DV) are based on a 2,000 calorie diet.

INGREDIENTS: Fractionated barley amylopectin (Vitargo® S2) [Source: Europe], citric acid, natural flavor, sucralose, natural color [beta carotene].

GLUTEN-FREE (via ELISA testing)

Manufactured in the USA at a cGMP facility.



Samples from each batch of Vitargo® S2™ multi-serving tubs are analyzed for substances banned by sport. For more details, see www.bscg.org

VitargoS2.com

877-GENR8-58 (877.436.7858)

Marketed exclusively by GENR8, Inc.

Dana Point, CA 92629

