

SUGGESTED USE: Shake lightly before each use. On an empty stomach or between meals, take 1 tsp (5 ml) by mouth. May mix in water or juice. Use as directed by your healthcare professional. Do not take within 3 hours of bedtime as it may interfere with sleep.

Children and pregnant or lactating women should consult their healthcare professional prior to use.

STORAGE: Store in a cool, dry place away from direct sunlight and out of reach by children. For maximum potency, after opening: refrigerate and use within 40 days.

Optimal Liposomal Vitamin C Plus provides 1,000 mg of Quali[®]-C vitamin C (sourced from Scotland) with rose hips extract and bioflavonoids naturally found alongside vitamin C in citrus fruits. The phospholipid base is non-hydrogenated and derived from GMO-free sunflower oil – without the common use of soy.

Quali[®]-C

Premium vitamins by DSM – www.qualityforlife.com
Quali[®] is a trademark of DSM

DOES NOT CONTAIN: Milk, eggs, fish, shellfish, tree nuts, peanuts, wheat, soy, GMO, artificial colors, or artificial flavors.

Optimal Liposomal Vitamin C Plus
Meets or Exceeds cGMP Quality Standards.



TARGETED SUPPORT[®]



Seeking Health[®]

Optimal Liposomal Vitamin C Plus

Dietary Supplement
5 FL. OZ. (150 ML)

Supplement Facts

Serving Size: 1 tsp (5 ml)

Servings Per Container: 30

	Amount Per Serving	%DV
Calories	10	
Total Carbohydrates	2 g	<1% [†]
Vitamin C	1,000 mg	1,111%
(as sodium ascorbate, ascorbic acid)		
Sodium	90 mg	4%
(as sodium ascorbate)		
Phosphatidylcholine and Phosphatidylethanolamine Complex (from non-GMO Sunflower oil)	450 mg	*
Bioflavonoid Blend	15 mg	*
Quercetin Dihydrate		
Rose Hips Extract		
Hesperidin		

*Daily Value (DV) not established.
[†]%DV based on a 2,000 cal/day diet.

Other Ingredients: Purified water, glycerin, natural flavors, gum arabic, and potassium sorbate (0.1%, to deter spoilage).

Find Other Effective Health Products At
www.seekinghealth.com

Manufactured for and Distributed by:
Seeking Health, LLC
3140 Mercer Ave
Bellingham, WA 98225
(800) 547-9812

VITACP-50Z-LIPO
REV. 010320

