

SUGGESTED USE: Take 2 capsules two times daily or as recommended by your health-care practitioner.

Tamper Evident: Use only if bottle is sealed. Store tightly sealed in a cool, dry place.

If pregnant, consult your health-care practitioner before using this product.


Manufactured for:

Tesseract Medical Research, LLC
1370 Trancas St., #376, Napa, CA 94558

www.TessMed.com

TetraCumin QR[®]

Curcumin Enhanced Delivery
Quick Release
Hypoallergenic



Tesseract
Medical Research[®]

Dietary Supplement
120 Vegetarian Capsules

SUPPLEMENT FACTS

Serving Size: Two Capsules

Servings Per Container: 60

Two Capsules Contain:	% DV
-----------------------	------

TetrahydroCurcumin Proprietary Complexes [†]	800 mg *
--	----------

*Daily Value (DV) not established.

Other Ingredients: Hypromellose (derived from cellulose) capsule, microcrystalline cellulose, calcium laurate, leucine, and plant derived enzyme.

[†]Solubility enhanced TetrahydroCurcumin (from *Curcuma longa* (rhizomes) in CyLoc[™], dextrin fiber matrix molecular traps. CyLoc[™] complexes compounded with DexKey[™] release agents (patent pending, AI Czap, Czap R&D, Licensed by Tesseract Medical Research, LLC).