

Directions: As a dietary supplement, mix 1 serving (2 level scoops) in 20 ounces of water and shake vigorously for 30 seconds. Sip throughout workout or athletic event.

An additional serving may be taken immediately after training or any time additional aminos are needed.

Using very cold water can hinder solubility and cause foaming.

Warning: This product is only intended to be consumed by healthy adults 18 years of age or older. Pregnant or nursing women should not use this product. Discontinue use and immediately consult your health care professional if you experience any adverse reaction to this product. Do not exceed recommended serving. Do not use if safety seal is broken or missing. **KEEP OUT OF REACH OF CHILDREN.**

† These statements have not been evaluated by the Food and Drug Administration. This product is not intended to treat, diagnose, cure or prevent any disease.

35 WORKOUTS

ULTRA MICRONIZED

Modern BCAA™

The Better BCAA™

MICRONIZED AMINOS

Are you still using dated, bulky BCAAs with non-micronized aminos & ordinary L-glutamine – loaded with unnecessary artificial flavors & dyes?

Get out of the stone age with USPlabs Modern BCAA™ – Equipped with 10 grams per serving of ultra-fine micronized aminos and zero artificial flavors or dyes.†



MIXES CLEAR*

*When mixed according to directions

Amino Acid Supplement



Net Wt. 12.3 oz (350g)

NEW!

UNFLAVORED

ZERO
ARTIFICIAL COLORS OR FLAVORS



Modern BCAAs
Mixes Clear*

Elapsed Time:
90 Seconds



Other BCAAs
Incomplete Dissolution
Loaded with Unnecessary
Flavors & Dyes

Supplement Facts

Serving Size 2 Scoops (10.00g)
Servings per Container: 35

	Amount per serving	% Daily Value
Proprietary Blend	10.00g	*
Micronized L-Leucine, Micronized L-Isoleucine, Micronized L-Valine, L-Alanyl-L-Glutamine (Sustamine™)		

* Daily value not established

SUSTAMINE™
Rehydrate. Replenish. Recover.™

Sustamine is a registered trademark of Kyowa Hakko Bio Co., LTD. Corporation

Ultra-fine & voluminous powder is a result of advanced micronization process.

Manufactured by USPlabs, LLC (2321 Marlene Dr., Ste 120, Dallas, Texas 75229) using Good Manufacturing Practices.

Questions and comments can be directed to our forum at www.usplabsdirect.com 1-800-850-3067



USA001