

PRODUCT BENEFITS: The Fish Oil supplied in PhysioLogics' EPA/DHA contains the polyunsaturated long-chain fatty acids Eicosapentaenoic Acid (EPA) and Docosahexaenoic Acid (DHA). Each serving delivers 300 mg EPA and 200 mg DHA to promote cardiovascular health.* Fish Oil supplies a concentrated source of these important fatty acids, which may be lacking from American diets.

Made from the oils of cold water, deep sea fish.
Marine Lipid Concentrate.

Directions: For adults, take one (1) softgel three times daily, preferably with a meal or follow the advice of your health care professional. As a reminder, discuss the supplements and medications you take with your health care providers.

WARNING: If you are pregnant, nursing, taking any medications or planning any medical procedure, consult your doctor before use. Discontinue use and consult your doctor if any adverse reactions occur.

*These statements have not been evaluated by the FDA. This product is not intended to diagnose, treat, cure or prevent any disease.

PROD. NO. 55217

To refill please call 1-800-765-6775
or visit www.physiologics.com

Carefully Manufactured by PhysioLogics, LLC
One Nutrition Plaza, Carbondale, IL 62901 U.S.A.



B51018

09G

PhysioLogics®

EPA | **DHA**
300 mg | 200 mg

Omega-3 600 mg

Promotes Cardiovascular Health*

Dietary Supplement
90 Softgels



Supplement Facts

Serving Size 1 Softgel

Amount Per Serving	%Daily Value	
Calories	10	
Calories from Fat	10	
Total Fat	1 g	2%**
Polyunsaturated Fat	<1 g	***
Cholesterol	10 mg	3%
Fish Oil	1,000 mg (1 g)	***
provides 600 mg of Omega-3 Fatty Acids, comprising of:		
EPA (Eicosapentaenoic Acid)	300 mg	***
DHA (Docosahexaenoic Acid)	200 mg	***
Other Omega-3 Fatty Acids	100 mg	***

**Percent Daily Values are based on a 2,000 calorie diet.

***Daily Value not established.

Other Ingredients: Gelatin, Glycerin, Mixed Natural Tocopherols.
Contains fish (anchovy, mackerel, sardine) ingredients.

FREE OF: artificial color, artificial flavor, artificial sweetener, preservatives, sugar, starch, milk, lactose, gluten, wheat, yeast, shellfish, sodium.

Purified to eliminate mercury

TAMPER RESISTANT: Do not use if seal under cap is broken or missing.
KEEP OUT OF REACH OF CHILDREN. Store in a dry place and avoid excessive heat.

Laboratory Tested to Guarantee Quality
For Health Care Professional Use Only
©PhysioLogics, 2011