
Suggested Use: Mix 1 to 2 scoops of Probolic-SR with 8 to 16 oz. of water in a shaker bottle or blender. Consume 2-3 servings a day or use as recommended by a physician or licensed nutritionist.

Sup	plem	ent	<b>Facts</b>
Serving	Size: 1 s	coop (33	g)
Serving	s Per Con	tainer:	•
Approx	imately 50	5	
Amount	Per Serving	9	

Servings Fer Container.		
Approximately 56		
Amount Per Serving		
Calories	140	
Calories from Fat	35	
%Daily Value		
<b>Total Fat</b> 4 g	6%†	
Saturated Fat 1.5 g	8%†	
<i>Trans</i> Fat 0 g		
Cholesterol < 5 mg	1%†	
<b>Total Carbohydrate</b> 5 g	2%†	
Sugars 0 g		
Dietary Fiber 2 g	8%†	
Protein 20 g	40%†	
Calcium 670 mg	67%	
Iron 4.3 mg	24%	
Phosphorus 520 mg	52%	
Sodium 310 mg	13%	
Potassium 510 ma	15%	

Amino Acid	<b>Profile</b>		
Mg Per 100g of Protein			
Alanine	4,230mg		
Arginine	7,040mg		
Asparatic Acid	11,130mg		
Cystine	1,250mg		
Glutamine*	21,710mg		
Glycine	3,830mg		
Histidine	2,600mg		
Isoleucine	5,950mg		
Leucine	7,650mg		
Lysine	6,500mg		
Methionine	1,380mg		
Phenylalanine	5,100mg		
Proline	5,430mg		
Serine	5,180mg		
Threonine	3,890mg		
Tryptophan	1,280mg		
Tyrosine	3,860mg		
Valine	6,010mg		

\*These statements have not been evaluated by the Food and Drug Administration. This product is not intended to diagnose, treat, cure or prevent any disease.





- Protect from heat, light and moisture

### Keep out of reach of children

- Store at 15°-30°C (59°-86°F).
- Do not purchase if seal is broken.





Anabolic/Anti-Catabolic Bio-Efficient Protein



# Sustained Release Micro-Feed Technology™

- Fast, Medium & Slow Release Profile Ideal Anytime: Pre-Workout, Post-Workout and Before Bedtime\*
- > Highest Levels of the Critical 5 Amino Acids (Glutamine, Arginine & BCAAs)\*
- Flax, CLA & MCTs

chocolate naturally and artificially flavored

## PROBOLIC-57—12 Hour Sustained Release Bio-Efficient Protein


Probolic-SR's patented 12 hour Sustained Release Micro-Feed Technology and superior Critical Flve Amino Acids Score (CFAAS) are major breakthroughs in protein supplementation. These two exculsive breakthroughs have set new standards in the bioefficiency and anabolic effects of protein by giving you something no other protein can: a continual 12 hour supply of the most critical amino acids to stimulate muscle growth, improve recovery—and prevent catabolic muscle wasting.\*

Amino acids are the building blocks for muscle tissue. In order to achieve maximum muscle growth, you must have an adequate supply of amino acids to feed your muscles. Certain amino acids are more important than others for stimulating and supporting muscle growth. These amino acids are called the Critical Five amino acids, which include glutamine, arginine and BCAAs (leucine isoleucine and valine). Probolic-SR's engineered protein matrix of whey concentrate, soy isolate and casein is precisely formulated to supply the abolute highest levels of these Critical Five amino acids. Aside from providing the highest Critical Five Amino Acid Score, this protein matrix also provides a fast, medium and slow "release profile", which improves bio-availability and nitrogen retention.\*

The exclusive patented Sustained Release Micro-Feed technology further enhances the "release profile" and provides a continual supply of Critical Five amino acids for up to 12 hours! This controlled sustained release of amino acids to your muscles increases bio-availability and nitrogen retention and creates the ideal anabolic/anti-catabolic state for maximum muscle growth and prevents muscle breakdown. Probolic-SR's Micro-Feed technology makes it the most bio-efficient protein available.\*

Feed Your Muscles with PROBOLIC-SR—the 12 Hour Muscle Feeder!





## PROBOLIC-SR's Release Profile:

- Fast Release Proteins provided by Whey—ideal for post-workout.
- Medium and Slow Release Proteins to feed muscle tissue between meals.
- Patented Micro-Feed Sustained Release Technology to provide a steady supply of amino acids for 12 hours.

Distributed by: Maximum Human Performance, Inc., 21 Dwight Place, Fairfield, New Jersey 07004 USA Toll Free: 1.888.783.8844, www.getMHP.com