Scan for product info

Read the entire label and follow the directions carefully prior to use.

DIRECTIONS: Take one (1) softgel daily with food, or as recommended by a healthcare practitioner.

Ubiquinone CoQ10 uses a patent-pending system designed to work with body heat and intestinal pH to enhance absorption of fat-soluble nutrients. When activated by ingestion, the Ubiquinone CoQ10 oil droplets are released to form a uniform and soluble emulsion with qastric fluids.

Refrigeration of this product may result in temporary discoloration (crystalization) in some softgels. However, such temporary crystallization is normall and does not affect product quality or performance, and will reverse with warmino.

Store tightly closed in a cool, dry place.

WARNINGS:

- KEEP OUT OF REACH OF CHILDREN.
- DO NOT EXCEED RECOMMENDED DOSE.
- Do not purchase if outer seal is broken or damaged.
 When using nutritional supplements, please consult

 When using nutritional supplements, please consult with your physician if you are undergoing treatment for a medical condition or if you are pregnant or lactating.

*These statements have not been evaluated by the Food and Drug Administration. This product is not intended to diagnose, treat, cure, or prevent any disease.



LifeExtension°

Super-Absorbable CoQ10

(Ubiquinone) • 100 mg with *d*-Limonene



Promotes Cardiovascular Health*

Dietary Supplement



100 Softgels

Supplement Facts

Serving Size 1 Softgel

Amount Per Serving

% Daily Value

Ubiquinone (as Kaneka Q10°)
Cold-Pressed Orange oil (peel)

100 mg

[std. to 90% *d*-Limonene (100 mg)]

**Daily Value not established.

Other ingredients: ethyl oleate, polysorbate 80, gelatin, glycerin, purified water, d-alpha tocopherol, rosemary extract.

Manufactured for:
Quality Supplements and Vitamins, Inc.
Ft. Lauderdale, Florida 33309
info@lifeextension.com • www.lef.org
To report a serious adverse event or obtain product

information, contact 1-866-280-2852.

Kaneka Q10° is a registered trademark of Kaneka Corporation.

Q00950G