

Hylavision has been formulated to provide nutritional supplement support for overall eye health.**

Hyaluronic Acid

*helps maintain healthy fluid balance in the eye***

Lutein/Zeaxanthin

*supports retinal and macular health***

Bilberry

*supports healthy eye circulation***

Vitamins A, B2, C, E, and Zinc

*powerful anti-oxidant properties***

For best results, use in conjunction with a healthy diet and exercise regimen.

Manufactured by:

hyalogic™

Questions?

610 NW Platte Valley Dr.
Riverside, MO 64150

toll-free: 866-318-8484

www.hyalogic.com



hyalogic™

hylavision

supports healthy vision**

*with Hyaluronic Acid
Lutein/Zeaxanthin, Bilberry,
and other Vitamins*

**Net contents 120 capsules
dietary supplement**

Supplement Facts

Serving Size
2 capsules

Servings Per Container
60

Amount Per Serving	% Daily Value †	
Vitamin A (as beta carotene)	6600 IU	132%
Vitamin C (as ascorbic acid)	400 mg	666%
Vitamin E (as D-alpha tocopherol)	150 IU	500%
Zinc (as zinc acetate)	60 mg	400%
Riboflavin (vitamin B2)	10 mg	588%
Manganese (as manganese acid chelate)	10 mg	500%
Copper (as copper gluconate)	2 mg	100%
Selenium (as selenium amino acid chelate)	40 mcg	57%
Bilberry Extract 4:1 (fruit)	40 mg	*
Hyaluronic Acid	5 mg	*
Lutein/Zeaxanthin Complex (maigold)	4 mg	*

* Daily Value not established.

† Daily Values based on a 2,000 calorie diet.

Other Ingredients: Vegetable Cellulose, Vegetable Magnesium Stearate, Silica.

FREE OF: Gluten, Wheat, Yeast, Sugar, Sodium, Artificial Flavor, Sweetener, Preservatives and Color

Directions: As a dietary supplement, consume two (2) capsules once a day; preferably with food.

Warning: Consult your physician prior to using this product if you are pregnant, nursing, taking medication or have a medical condition. Tamper resistant. Do not use if seal is broken. Keep out of reach of children. Store in cool, dry place.

****These statements have not been evaluated by the Food and Drug Administration. This product is not intended to diagnose, treat, cure, or prevent any disease.**

HLVSN1113