V1.0



Made in the UK for Seroyal USA, Pittsburgh, PA 15275, 1-888-737-6925

Recommended intake: (Children 1-4 years): In a glass, add water or milk to one scoop of HLC Baby B and mix. Take once daily with a meal, at least 2-3 hours before or after taking antibiotics, or as professionally directed

For risk information, see outer box.

Vacuum-sealed for your protection and for product freshness. Do not use if outer seal is missing or broken. Store under refrigeration (below 46 °F). Keep out of reach of children

HLC BABY B

probiotic supplement for breast fed babies

Pack size: Net Wt 0.2 oz (6 g)

Supplement Facts

Serving Size 200 mg, Serving 30, Amount per serving: HLC Consortium (10 billion viable cells t) of: Bifidobacterium bifidum (CUL-20); Bifidobacterium animalis subsp. lactis (CUL-34): Lactobacillus salivarius (CUL-61); Lactobacillus paracasei (CUL-08) † % Daily Value (DV) not established

Supplement Facts

Bifidobacterium bifidum (CUL-20)

animalis subsp. lactis (CUL-34)

10 billion viable cells t

Serving Size 1 Scoop (200 mg)

Servings per Container 30

Each Serving Contains

HLC Consortium

Bifidobacterium

Contains: Milk

Recommended intake: Children (1-4 years): In a glass, add water or milk to one scoop of HLC Baby B and mix. Take once daily with a meal, at least 2-3 hours before or after taking antibiotics, or as professionally directed.

Risk information: If you have an immune-compromised condition. do not use. If you are experiencing fever, vomiting, bloody diarrhea or severe abdominal pain, consult your healthcare practitioner prior to use. If symptoms of digestive upset occur. worsen or persist beyond three days. discontinue use and consult your healthcare practitioner.

Safety-sealed for your protection and for product freshness. Store under refrigeration (below 46 °F).

Keep out of reach of children.

HLC BABY B probiotic supplement for breast fed babies



Human Lactic Commensals (HLC) probiotic strains were studied in a randomized. double-blind, placebo-controlled human trial at 10 billion CFU per day, Proprietary non-pathogenic human-sourced microflora bacterial strains are utilized as they have strong epithelial adherence and a naturally high tolerance to stomach acid.

Lactobacillus salivarius (CUL-61) Lactobacillus paracasei (CUL-08) † % Daily Value not established

Other ingredients: Transgalactooligosaccharides (GOS) (from milk)

BARCODE LOCATION 883196 21670

Code PR28-6

Made in the UK for Seroval USA Pittsburgh, PA 15275, 1-888-737-6925

Pack size: Net Wt 0.2 oz (6 g)