incentrated Series

The Super Concentrated Series are the most powerful to constraint of the Concentration of the Concentrated Series in precise and fully disclosed clinically research of the Concentrated Series is powered by premium quality ingredients in procise and fully disclosed clinically researched doses. No propriets when the Concentrated Series is powered to the Concentrated Series in precise and fully disclosed clinically researched doses. No procise and control of the Concentration of the Concent

The Most Powerful Nitric Oxide Stimulant Ever.

NanoStim is a powerful super-concentrated nitric oxide stimulant hat delivers superior muscle pumps for an incredible training experience. Based on the latest clinical research on nitric oxide metabolism, NanoStim is formulated with the scientifically studied compounds L-citrulline, quercetin, yohimbe extract and huperzine A.A.

Why You Should Use NanoStim¹¹

 Unlike many other nitric oxide products, NanoStim™ contains a powerful 3000mg dose of L-otirulline, which is shown in a study on nitric oxide metabolism to be superior to a larger dose of arginine for increasing blood levels of arginine.

Schwedhelm et al., 2007. Br Journal of Clinical Pharmacology, 65(1):51-59

NanoStim™ delivers an extremely potent dose of yohimbe extract, supplying both rauwolscine and yohimbine. These have been shown to block alpha-2 receptors, which can play a significant role in regulating blood flow.

Formulated with quercetin, this powerful compound has been clinically shown to increase the key markers of nitric oxide

production and reduce a vasoconstricting factor.▲ Loke et al., 2008. American Journal of Clinical Nutrition. 88:1018-1025.

Huperzine-A has the ability to block the enzymes that decrease acetylcholine levels. Increasing acetylcholine levels can improve focus.

100% Money-Rack Guarante

If this is not the most powerful nitric oxide stimulant you've ever tried, simply return the unused portion with your receipt to place of





NANOSTIM

Super-Concentrated Nitric Oxide Stimulant*

Pump • Focus • Vascularity*



(as Huperzia serrata) (aerial parts + Daily Value not established

OTHER INGREDIENTS DICA CIUM PHOSPHATE HYDROXYPROPYLCELLILOSE, COATING (HYPROMELLOSE, TITANUM, DIOXIDE, POLYDEXTROSE, TALC, MALTODEXTRINDEXTRIN, MEDIUM, CHAN, TRIGLYCERDES,

100mca

Directions: Take 2 caplets of NanoStim™ with a glass of water 30-45 minutes before your workout to assess tolerance. For full effects, take 4 caplets before your workout. Do not exceed 4 caplets in a 24-hour period. Read the entire label before use and follow directions provided.

WARNING: Not intended for use by persons under 18. Do not use if pregnant or runsing, Discontinue size and consult a medical doctor if you experience unusual symptoms, Consult a medical doctor before use if you have been frested for or diagnosed with, or lave a laminy instruy of you have been frested for or diagnosed with, or lave a laminy instruy of or heat problems or if you are using any prescription or over-the-counter drugis, including blood thinners. Do not combine with other stimulants. Consult a medical doctor before starting any det or exercise program. Do not exceed recommended serving, improper use of the product, will not exceed recommended serving, improper use of the product, will not exceed recommended serving, improper use of the product, will not exceed recommended serving, improper use of the a count, or the control of the control of the control of the control of the counter of the counter of the control of the control of the control of the counter of the counter

Patent pending, Made in the U.S.A. from domestic and international ingredients. Distributed by lovate Health Sciences U.S.A. Inc. 1105 North Market Street, Suite 1330, Wilmington, DE 19801. © 2011. For lot no. and expiry date: see bottle.

EOR POSITION ONLY 31656 6030111

cc